

UltraPan 100 Developing times

Developer	Dilution	Time	Temperature
D76/Flic Classic MQ	Stock	6-7	20c
D76/Flic Classic MQ	1+1	10	20c
Flic B/W & Green	Stock	11.75	20c
ID-11	Stock	6-7	20c
Ifosol 3	1+9	5	20c
Ifosol 3	1+14	7.5	20c
HC-110	B	6	20c
HC-110	F	12	20c
HC-110	H	9	20c
Rodinal	1+50	8-9	20c
Rodinal	1+100	18	20c
Xtol	Stock	5-6	20c
Xtol	1+1	8	20c

UltraPan 400 Developing times

Developer	Dilution	Time	Temperature
D76/Flic Classic MQ	Stock	7-8	20c
D76/Flic Classic MQ	1+1	11-12	20c
Flic B/W & Green	Stock	11.75	20c
ID-11	Stock	7-8	20c
ID-11	1+1	11-12	
Ifosol 3	1+9	8.5	20c
HC-110	B	7	20c
HC-110	H	13	20c
Rodinal	1+50	11.75	20c
Rodinal	1+100	30	20c
Xtol	Stock	7	20c
Xtol	1+1	9.5	20c